


OsiriX Radiology Solutions with Apple



Ultrafast performance



Intuitive interactive user interface



FDA & CE certified

OsiriX MD is the complete edition of OsiriX. With ultrafast performance and an intuitive interactive user interface, it is the **most widely used medical images viewer in the world**. It fully supports the DICOM standard for an easy integration in your workflow environment and an open platform for development of processing tools. It offers **advanced post-processing techniques in 2D and 3D**, exclusive innovative techniques for **3D & 4D navigation**, including PET-CT and SPECT-CT fusion, and a complete integration with any PACS server. It can import and display any medical images (DICOM files) from CD, DVD, USB stick and others, OsiriX MD **supports 64-bit computing and multithreading** for the best performances on the most modern computers.

OsiriX MD is **FDA, CE and ANVISA** certified and validated for clinical use in medicine.

VisionMed OsiriX workstations have been designed keeping in mind the financial savings for radiology processing costs when using **Apple** computer-based systems. It can help small hospitals to get into PACS with minimal cost. A small hospital that can't afford a PACS can access images through Mac, download OsiriX MD, and look at images. The powerful medical imaging software OsiriX MD with applications such as VMware Fusion and others can run any windows applications needed and is less expensive and better than any diagnostic workstation.

“The unique work flow management tools make OsirixPacs a truly outstanding PACS. Ease of use, reliability and remarkable pricing make OsirixPacs an obvious choice for hospitals and imaging centers.”

Andrew Lee, M.D., Chino Hills, CA

“A great PACS system with outstanding support. Osirix has successfully met the IT needs of our imaging facility.”

Sheena Shimma, Owner, Northwest Prenatal Ultrasound



INSPIRE MEDITECH PRIVATE LIMITED

407 Red Rose Building 49-50, Nehru Place New Delhi-110019
 Office +91 11 26284578, 79, 80
 H/P+919811078156, 919313899647, 919717480657
 sales@inspiremeditech.com | www.inspiremeditech.com