



VisionMed

HIGH DEFINITION RECORDER FOR SURGICAL PROCEDURES

Record Video the Easy Way!

VisionMed makes recording video easier than ever. You can record high definition video from virtually any medical video device including endoscopes, arthroscopes, surgical cameras, ultrasounds, etc. Records upto 1080p video into USB stick or Hard Disk by one click, PC not needed. With microphone input, it allows surgeon to record live voice from microphone with video together.

Highest Quality Video. VisionMed high performance hardware and advanced compression algorithms provide the medical industry's clearest videos and images. Every pixel of every frame is captured with the truest possible fidelity.

Universal Compatibility. VisionMed recorders can be paired with virtually all makes and models of medical cameras and monitors. Videos recorded with VisionMed can be viewed on any modern computer, laptop or tablet.

Compact File Size. The high-definition VisionMed also saves video in the H.264/MP4 format that offers the industry's highest quality HD images with the most compact file size.

Rugged Case. The unit is protected by a rugged metal case for use in demanding medical environments.

- Interface : USB 2.0 Host.
- Video Input : HDMI, Component Video, Composite.
- Storage Media : USB Flash Disk, Portable Hard Disk,
- Resolution : 1920 x 1080 (60p).
- Recording : H.264/MP4, schedule recording supported.
- Power Supply : USB B, DC 5V / 2A.



By combining the latest technology with ease-of-use and an affordable price tag, VisionMed recorders have become a "must have" tool for doctors and OEMs.

INSPIRE INFOTECH PRIVATE LIMITED
MEDICAL DIVISION

407 | Red Rose Building | 49-50 Nehru Place | New Delhi - 110019

Tel: +91 26284578-79 | 9811078156

medical@inspireinfo.com | sales@inspirehealthcare.in | www.inspirehealthcare.in

inspire
HEALTHCARE
THE DR PEOPLE